



DATAMARS

COMPACT MAX Universal Scanner

PetLink's Universal Scanner, the COMPACT MAX, features top of the line readings distances due to its peripheral antenna. The scanner can easily fit in a pocket and is able to read FDX-B (ISO), FDX-A/Fecava, Trovan, and Avid encrypted microchips. It also has basic water protection, making it an ideal scanner for any condition.

From all species, ages, and sizes, the COMPACT MAX delivers the identification information you need. The COMPACT MAX has great ergonomics and an intuitive, icon-based, graphic interface - which makes it the easier reader to use on the market.

Highlights

- Best-in-class reading performances
- Long-lasting, rechargeable lithium battery
- Slim and elegant design
- Extremely easy to use
- Bright color display
- Ideal for vet clinics, shelters, rescues, and animal control agencies

Keys to Effective Scanning:

- Hold the scanner close to or touching the animal.
- Scan slowly and repeatedly over the entire body to ensure a microchip is present.
- Scan first in a horizontal, then a vertical "S" pattern down the animal's body.
- Consider scanning each animal more than once to ensure any present microchip can be localized by the scanner.
- Remove metal collars before scanning - computers, metal tables, and even fluorescent lighting can interfere with effective scanning

Technical Specifications

Dimensions

- 6.29" x 3.14" x 0.78"

Weight

- 6 oz

Operating Temperature

- +23 F to +131 F 95% RH, non-condensing

Charging Temperature

- +41 F to +113 F 95% RH, non-condensing

Storage Temperature

- -4 F to +149 F 95% RH, non-condensing

Battery

- Li-ion (1250 mAh)

Typical Reading Distance

- FDX-B (134.2 kHz) 4.72"
- FDX-A/Fecava (125 kHz) 2.76"
- Trovan (128 kHz) 1.55"
- Avid Encrypted (125 kHz) 2.76"

Housing

- White ABS, water-resistant with USB cap inserted

Approved by:

- CE/FCC
- National approval in Canada, Australia, and New Zealand

Standards

- ISO 11784/5 (FDX-B) non-synchronizing transceiver

Frequency

- 134.2 kHz, 125 kHz, and 128 kHz

Questions?

Contact us at (781) 281-2216
or at ussales@petlink.net

