

| Brand | | ASC Recommended Screw Torque Value |
|---------------------|-------------------------------------|---------------------------------------|
| Dentis | OneQ-SL/S-Clean Mini | 15N. Cm |
| | OneQ-SL/S-Clean Regular | 25N. Cm |
| | i-clean Regular | 25N. Cm |
| Dentium | SuperLine (Implantium) Regular | 25N. Cm |
| | SimpleLine II Regular | 25N. Cm |
| | SimpleLine II Wide | 25N. Cm |
| Dentsply Sirona | Astra Tech™ OsseoSpeed® EV 3.0 | 25N. Cm |
| | Astra Tech™ OsseoSpeed® EV 3.6 | 25N. Cm |
| | Astra Tech™ OsseoSpeed® EV 4.2 | 25N. Cm |
| | Astra Tech™ OsseoSpeed® EV 4.8 | 25N. Cm |
| | Astra Tech™ OsseoSpeed® EV 5.4 | 25N. Cm |
| | Astra Tech™ OsseoSpeed® TX 3.0 | 15N. Cm |
| | Astra Tech™ OsseoSpeed® TX 3.5, 4.0 | 20N. Cm |
| | Astra Tech™ OsseoSpeed® TX 4.5, 5.0 | 25N. Cm |
| | Ankylos C/X | |
| DIO | UF Mini | 15N. Cm |
| | UF Regular | 25N. Cm |
| | UF Internal Regular | 25N. Cm |
| | UF Internal Wide | 25N. Cm |
| | SM Small | 25N. Cm |
| | SM Regular | 25N. Cm |
| Hiossen (Osstem) | Hiossen ET (Osstem ET) Mini | 20N. Cm |
| | Hiossen ET (Osstem ET) Regular | 25N. Cm |
| | Hiossen HG Ultra-Wide | 25N. Cm |
| | Hiossen HS (Osstem SS) Regular | 25N. Cm |
| | Hiossen HS (Osstem SS) Wide | 25N. Cm |
| Megagen | AnyOne® Mini | 15N. Cm |
| | AnyOne® Regular | 25N. Cm |
| | ExFeel | 25N. Cm |
| MIS | M4/Seven Standard | 25N. Cm |
| | M4/Seven Wide | 25N. Cm |
| Neobiotech | IS Regular | 25N. Cm |
| | IT Regular | 25N. Cm |
| | IT Wide | 25N. Cm |
| Nobel Biocare | NobelActive™ NP | 25N. Cm |
| | NobelActive™ RP | 25N. Cm |
| | NobelActive™ WP | 25N. Cm |
| | NobelReplace™ NP Conical Connection | 25N. Cm |
| | NobelReplace™ RP Conical Connection | 25N. Cm |

| Manufacturer | Brand | ASC Recommended Screw Torque Value |
|-------------------|--|------------------------------------|
| Straumann® | Bone Level 3.3 NC (Narrow CrossFit™) | 25N. Cm |
| | Bone Level 4.1, 4.8 RC (Regular CrossFit™) | 25N. Cm |
| | synOcta 4.8 RN | 25N. Cm |
| | synOcta 4.8 WN | 25N. Cm |
| Zimmer® Dental | Zimmer Eztetic™ 3.1 | 25N. Cm |
| | Tapered Screw - Vent™ 3.5, 4.5, 5.7 | 25N. Cm |
| | Biomet 3i Certain™ 3.4 | 20N. Cm |
| | Biomet 3i Certain™ 4.1 | 20N. Cm |
| | Biomet 3i Certain™ 5.0 | 20N. Cm |
| | Biomet 3i Certain™ 6.0 | 20N. Cm |

| Manufacturer | Brand | ASC Recommended Screw Torque Value |
|-----------------|--|------------------------------------|
| Implant Direct | Legacy 3.0 | 25N. Cm |
| Keystone Dental | PrimaConnex™/Genesis™ 3.8 SD | 25N. Cm |
| | PrimaConnex™/Genesis™ 4.5 RD | 25N. Cm |
| | PrimaConnex™/Genesis™ 5.5/6.5 WD | 25N. Cm |
| Neoss | ProActive Straight 3.25 | 25N. Cm |
| | ProActive Straight 3.5, 4.0, 4.5, 5.0, 5.5 | 25N. Cm |
| | ProActive Tapered 3.5, 4.0, 4.5, 5.0, 5.5 | 25N. Cm |
| Nobel Biocare | NobelActive™ 3.0 | 15N. Cm |
| Straumann® | Bone Level 2.9 SC (Narrow CrossFit™) | 25N. Cm |