



DESIGNED IN
COLLABORATION WITH:



AARON DALUISKI, MD
HOSPITAL FOR SPECIAL SURGERY, NEW YORK, NY

5613 CAST-AWAY®



INDICATIONS:

- STABLE ELBOW FRACTURES
- ELBOW DISLOCATION
- DISTAL BICEPS REPAIR AND OTHER SURGICAL INTERVENTIONS
- ELBOW INJURY

FEATURES AND BENEFITS:

- **DUAL ROD STAYS** – PROVIDE SUPERIOR RIGIDITY AND SUPPORT WHILE REMAINING LOW PROFILE
- **LIGHT-WEIGHT, LOW-PROFILE DESIGN** – LEAVES THE ADJACENT, UNAFFECTED JOINTS UNRESTRICTED AND PROMOTES COMPLIANCE IMPROVING CLINICAL OUTCOMES
- **SIMPLE APPLICATION** – PROMOTES GOOD HYGIENE/SKIN CARE AND ELIMINATES CAST COMPLICATIONS
- **VERSATILE STRAPPING SYSTEM** – ALLOWS FOR CUSTOMIZED FIT AND COMFORT
- **TELESCOPING HAND CUFF** – SUPPORTS THE WRIST AND CONTROLS FOREARM ROTATION WHEN NEEDED
- **REMOVABLE SLING STRAP** – ELIMINATES THE EARLY NEED FOR A SLING IN CONJUNCTION WITH THE BRACE
- **ADJUSTABLE HEIGHT HUMERAL CUFF** – ALLOWS FOR A CUSTOMIZED FIT



OPTIONAL ACCESSORY:
ORTHOWEAR ELBOW (5614)
* SIZES XXS THROUGH XL

USPN 9,782,285



5613 ORDERING INFORMATION:

SIZE	ARM CIRCUMFERENCE (IN)	H (IN)
S	8-10	8
M	10-12	9
L	12-14	10

MEASURE WIDEST PART OF FOREARM
HEIGHT OF HUMERUS INCLUDED FOR REFERENCE

CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM



HELY & WEBER
BRACING INNOVATIONS



© 2024 HELY & WEBER