### Hearing aid selection guide.

## **TruHearing**®

#### WHAT MATTERS MOST TO YOU?

Choose all that apply

	Hearing in challenging				
environments					



Hear conversations more easily with greater focus on speech and less on the challenging soundscape around you.





Hear your own voice as nature intended - avoid having your voice sound too loud or unnatural. Increase your acceptance and keep engaged in conversations.

#### Connectivity



Stream calls, TV, video, music, and more to both ears. Easy to use - works with TruHearing App.

#### Rechargeability



Perfect for those who want to avoid changing small hearing aid batteries, especially if you have dexterity issues.

#### Tinnitus relief



Therapy strategies provide a soothing, positive listening experience.

#### Personalized sound



Hear sounds and speech from any direction - automatically adjusts to adapt to your unique listening situations and movements.

#### **STYLE OPTIONS**









RECEIVER IN CANAL (RIC)



BEHIND THE EAR (BTE)



READY-TO-WEAR



CUSTOM IN THE EAR (ITE)

#### **INCLUDED WITH PRODUCT**



The all-in-one app works not only as a remote control for your hearing aids, but a control of Bluetooth streaming, and even helps you track your fitness, steps, and hearing aid use.



#### **NOTES**

#### PROFESSIONAL RECOMMENDATIONS

	Select which is best for your:		
Hearing aid options	Hearing	Style	

# Patient Benefits by product and performance level.



Performance Level	Premium	Advanced
Hear conversations more easily with greater focus on speech and less on the challenging soundscape around you. (RIC, BTE)	J	J
Automatic adjustments via Motion Sensors to adapt to unique listening needs the environment changes. (RIC, BTE)	J	J
Natural voice and higher acceptance via Own Voice Processing to avoid having your voice sound too loud or unnatural. (RIC, BTE)	J	J
Hearing phone calls in both ears by using the TwinPhone feature- holding up phone to ear and having it stream automatically to the opposite ear OR take advantage of Bluetooth® (BT) streaming via a smartphone to hear the phone conversation in both ears*. (RIC, BTE, Instant-fit, Custom)	J	J
Enhance music with preset programs via HD Music. Included options for live music, recorded music, and personal performances. (RIC, BTE, Instant-fit, Custom)	3	3
Relief for ringing in the ears (Tinnitus) via Therapy solutions and exclusive Notch Therapy. (RIC, BTE, Custom)	J	J
Automatically focus on speakers in front of you via Directionality for better speech understanding in noise. (RIC, BTE, Instant-fit, Custom)	J	J
Customize focus in the TruHearing App via Directional Hearing to adjust as needed to better suit your personal preferences in any situation - dining in restaurants, walking outdoors, or visiting with friends. (RIC, BTE, Custom)	Front Back Right Left	Front Back Right Left
Automatic spatial awareness via Spatial SpeechFocus to hear speech and noise from all directions, like hearing those around you when in a car. (RIC, BTE, Custom)	Front Back Right Left	Front Back Surround
Reduced wind noise via eWindscreen for better clarity in speech when outdoors on a windy day or when riding in a car with the top down. (RIC, BTE, Instant-fit, Custom)	Auto	Manual
Better awareness of high-frequency sound via Extended Bandwidth to hear higher pitches with enhanced sound quality, such as a bird chirping, a small child or an alarm. (RIC, BTE, Instant-fit, Custom)	J	
Increase understanding in reverberant environments via echoShield to reduce echoes in large spaces like places of worship, museums, and auditoriums. (RIC, BTE, Instant-fit, Custom)	J	

Higher performance levels include more channels for additional customization options. Checkmarks highlight feature availability at respective levels, while numbers / text highlight added functionalities.