

640 MTC FRACTURE BRACE™

A unique alternative to a traditional sugar tong splint or Munster Cast.



FEATURES AND BENEFITS:

- **TWO-PIECE COMPONENT DESIGN:**
 - Allows the forearm to be placed in any degree of supination or pronation for optimal healing
 - Makes fitting a variety of arm lengths easy
 - Lends itself to a simple step-down (wrist only) application as healing progresses
- **DISTAL COMPONENT:**
 - has volar and dorsal stays providing excellent wrist stabilization
- **OPTIONAL D-RING STRAP:**
 - a versatile option to ensure a secure fit
- **LIGHTWEIGHT, BREATHABLE MATERIALS:**
 - make the brace especially comfortable compared to heavy casting or ridged splinting

INDICATIONS:

- **STABLE WRIST FRACTURES**
- **TFCC INJURY/REPAIR**
- **ULNAR OSTEOTOMY PROCEDURES**
- **MULTI-TRAUMA INJURIES**
- **POST-OPERATIVE IMMOBILIZATION**



**EFFECTIVE.
CONVENIENT.
COMFORTABLE.**

DESIGNED IN
COLLABORATION WITH:



JAMIE MCMILLAN, OTR, CHT



**PRIMARY GOAL:
PREVENT FOREARM
ROTATION WHILE ALLOWING
FUNCTIONAL ELBOW RANGE
OF MOTION.**



DISTAL COMPONENT DIMENSIONS:

640
XS (18.4cm)
8" (20cm)
10" (25cm)



PROXIMAL COMPONENT DIMENSIONS:

	SMALL	LARGE
W	9.7" (25cm)	12" (30cm)
L	8" (20cm)	9" (23cm)



640 ORDERING INFORMATION:

SIZE	ARM LENGTH	ARM CIRCUMFERENCE	COMPONENTS USED
640-S-8-XS	10"-14" (25-36cm)	7.5"-11" (19-28cm)	8" XS WRIST: S CUFF
640-S-8	10"-15" (25-38cm)	9"-11" (23-28cm)	8" WRIST: S CUFF
640-S-10	11"-17" (28-43cm)	9"-11" (23-28cm)	10" WRIST: S CUFF
640-L-8	11"-16" (28-41cm)	11"-14" (28-36cm)	8" WRIST: L CUFF
640-L-10	11"-18" (28-46cm)	11"-16" (28-41cm)	10" WRIST: L CUFF

BOLD DENOTES MOST POPULAR SIZES

SPECIFY RIGHT OR LEFT

ARM LENGTH: MEASURE FROM OLECRANON TO DISTAL PALMAR CREESE

ARM CIRCUMFERENCE: MEASURE WIDEST PART OF FOREARM

CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM



HELY & WEBER
BRACING INNOVATIONS



© 2026 HELEY & WEBER